

Anti-Ferritin heavy chain 1/FTH1 Antibody (7Y790)

Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	7Y790
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:1000-1:2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human Ferritin heavy chain 1 / FTH1 Protein (TMPY-02711)
Antigen Species:	Human
Synonyms:	PLIF;HFE5;ferritin, heavy polypeptide 1;FTHL6;FHC;FTH;PIG15

Research Background

FTH1 (ferritin, heavy polypeptide 1) is the heavy subunit of ferritin which is the major intracellular iron storage protein in prokaryotes and eukaryotes. It is composed of 24 subunits of the heavy and light ferritin chains. Variation in ferritin subunit composition may affect the rates of iron uptake and release in different tissues. A major function of ferritin is the storage of iron in a soluble and nontoxic state. Defects in ferritin proteins are associated with several neurodegenerative diseases. FTH1 gene has multiple pseudogenes. Several alternatively spliced transcript variants have been observed, but their biological validity has not been determined. FTH1 stores iron in a soluble, non-toxic, readily available form. It is important for iron homeostasis. It has ferroxidase activity. Iron is taken up in the ferrous form and deposited as ferric hydroxides after oxidation. It also plays a role in delivery of iron to cells. FTH1 mediates iron uptake in capsule cells of the developing kidney.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481