

Anti-p63 Antibody (3I778)

Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	3I778
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:1000-1:2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human p63 / TP63 protein (TMPY-02651)
Antigen Species:	Human
Synonyms:	LMS;p53CP;OFC8;p40;RHS;TP53CP;p73H;EEC3;SHFM4;B(p51B);p63;TP73L;NBP;KET;p51;tumor protein p63;p73L;TP53L;AIS;B(p51A)

Research Background

Tumor protein p63 is a protein also known as transformation-related protein 63, TP63, and p63. Tumor protein p63 / p63 is a member of the p53 family of transcription factors whose members P53, p63, and p73 have similar features in their gene structures and functions. An animal model, p63^{-/-} mice has been useful in defining the role p63 plays in the development and maintenance of stratified epithelial tissues. This p63 encoding protein p63 has a dramatic impact on replenishment of cutaneous epithelial stem cells and on ovarian germ cell survival. Although these two fundamental roles of p63 attest to its powerful place in development, its other functions, specifically the apparent capacity of p63, is to supervise the emergence of new cell populations in the breast, prostate, cervix, and upper reproductive tract. P63^{-/-} mice have several development defects which include the lack of limbs and other tissues, such as teeth and mammary glands, which develop as a result of interactions between mesenchyme and epithelium. Mutations in this protein are associated with ectodermal dysplasia, and cleft lip / palate syndrome 3, ADULT syndrome (acro-dermato-ungual-lacrimal-tooth), limb-mammary syndrome, et al.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481