

Anti-BPIFB1 Antibody (5H938)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	5H938
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:5000-1:10000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human BPIFB1 / LPLUNC1 Protein (TMPY-03851)
Antigen Species:	Human
Synonyms:	C20orf114;BPI fold containing family B, member 1;LPLUNC1

Research Background

BPIFB1, also known as LPLUNC1, belongs to the BPI/LBP/Plunc superfamily, plunc family. BPIFB1 may be involved in the innate immune response to bacterial exposure in the mouth, nasal cavities, and lungs. BPIFB1 is expressed in the upper respiratory tract and oral cavity, which may function in host defence. The expression of BPIF proteins is associated with CF lung disease in humans and mice. It is unclear if this elevation of protein production, which results from phenotypic alteration of the cells within the diseased epithelium, plays a role in the pathogenesis of the disease. BPIFB1 is an abundant, secreted product of goblet cells and minor mucosal glands of the respiratory tract and oral cavity and suggest that the protein functions in the complex milieu that protects the mucosal surfaces in these locations.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481