

## Anti-Complement C7 Antibody (2B681)

### Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	2B681
Purification:	Protein A

### Applications

Application:	ELISA(Cap)
--------------	------------

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human Complement component 7 Protein (TMPY-02892)
Antigen Species:	Human
Synonyms:	Complement C7; complement component 7

### Research Background

Complement component 7 is a component of the complement system. It belongs to the complement C6/C7/C8/C9 family. It contains 1 EGF-like domain, 1 LDL-receptor class A domain, 1 MACPF domain, 2 Sushi (CCP/SCR) domains and 2 TSP type-1 domains. Complement component 7 serves as a membrane anchor. It participates in the formation of Membrane Attack Complex (MAC). People with C7 deficiency are prone to bacterial infection. It is a constituent of MAC that plays a key role in the innate and adaptive immune response by forming pores in the plasma membrane of target cells. Defects in C7 are a cause of complement component 7 deficiency (C7D). A rare defect of the complement classical pathway associated with susceptibility to severe recurrent infections, predominantly by *Neisseria gonorrhoeae* or *Neisseria meningitidis*.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481