

Anti-Vasorin Antibody (5N124)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	5N124
Purification:	Protein A

Applications

Application:	ELISA(Cap)
--------------	------------

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human Vasorin Protein (TMPY-03058)
Antigen Species:	Human
Synonyms:	Protein Slit-Like 2;VASN;Vasorin;SLITL2

Research Background

Vasorin is a typical type I membrane protein. It contains 1 EGF-like domain, 1 fibronectin type-III domain, 10 LRR (leucine-rich) repeats, 1 LRRCT domain and 1 LRRNT domain. Vasorin is predominantly expressed in vascular smooth muscle cells, and that its expression is developmentally regulated. Vasorin directly binds to transforming growth factor (TGF)- β and attenuates TGF- β signaling in vitro. This suggests that down-regulation of vasorin expression contributes to neointimal formation after vascular injury and that vasorin modulates cellular responses to pathological stimuli in the vessel wall.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481