

## Anti-CFHR2 Antibody (9T861)

## Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	9T861
Purification:	Protein A

## Applications

Application:	ELISA
Recommended	ELISA: 1:5000-1:10000

## Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	Recombinant Protein: Human CFHR2 Protein (TMPY-03061)
Antigen Species:	Human
Synonyms:	FHR2; complement factor H-related 2; CFHL2; HFL3
Biology Area:	Serine Proteases and Regulators

## Research Background

CFHR2 belongs to the complement factor H protein family. The human complement factor H protein family consists of the complement and immune regulators factor H, the factor H-like protein 1 (FHL-1) and five factor H-related proteins (CFHR-1 to -5). Members of the H-related protein family are exclusively composed of individually folded protein domains, termed short consensus repeats (SCRs) or complement control modules. CFHR2 contains 4 Sushi (CCP/SCR) domains and is expressed by the liver and secreted in plasma. CFHR2 might be involved in complement regulation. It can associate with lipoproteins and may play a role in lipid metabolism.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481