

## Anti-CD177 Antibody-PE (3S149)

### Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	PE
Clone:	3S149
Purification:	Protein A

### Applications

Verified Activity:	Flow cytometric analysis of Human CD177 expression on human whole blood granulocytes. Cells were stained with PE-conjugated anti-Human CD177. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable granulocytes.
Application:	FCM
Recommended	5 µl/Test, 0.1 mg/ml

### Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight. Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human CD177 / PRV-1 / PRV1 Protein
Antigen Species:	Human
Synonyms:	NB1GP;HNA2A;PRV-1;NB1;PRV1;HNA-2a;CD177 molecule

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481