

## Anti-TNFAIP8 Antibody (30782)

## Product Details

Ig Type:	Mouse IgG2a
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	30782
Purification:	Protein A

## Applications

Application:	ELISA
Recommended	ELISA: 1:1000-1:2000

## Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	Recombinant Protein: Human TNFAIP8 Protein (TMPY-03416)
Antigen Species:	Human
Synonyms:	tumor necrosis factor, alpha-induced protein 8;MDC-3.13;SCCS2;NDED;SCC-S2;tumor necrosis factor, $\alpha$ -induced protein 8;GG2-1

## Research Background

Tumor necrosis factor (TNF)-alpha-induced protein 8 (TNFAIP8) family is a newly identified protein with vital roles in maintaining immune homeostasis. Tumor necrosis factor-alpha-inducible protein 8 (TNFAIP8) is a TNF-alpha inducible anti-apoptotic protein with multiple roles in tumor growth and survival. by the creation of cellular autophagy events, TNFAIP8 promotes cell survival and drug resistance in prostate cancer cells. TNFAIP8 regulates Hippo pathway through interacting with LATS1 to promote cell proliferation and invasion in lung cancer. TNFAIP8 may serve as a candidate biomarker for poor prognosis and a target for new therapies.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481