

Anti-DCUN1D2 Antibody (1X996)

Product Details

| | |
|---------------|--------------|
| Ig Type: | Rabbit IgG |
| Reactivity: | Human |
| Conjugation: | Unconjugated |
| Clone: | 1X996 |
| Purification: | Protein A |

Applications

| | |
|--------------------|---|
| Verified Activity: | <p>1. Anti-DCUN1D2 rabbit monoclonal antibody at 1:500 dilution.</p> <ul style="list-style-type: none">-Lane A: HCT116 Whole Cell lysate.-Lysates/proteins at 30 µg per lane.-Secondary-Goat Anti-Rabbit IgG H&L (Dylight800) at 1/10000 dilution.-Developed using the Odyssey technique.-Performed under reducing conditions.-Predicted band size:30 kDa.-Observed band size:30 kDa. <p>2. DCUN1D2 was immunoprecipitated using:</p> <ul style="list-style-type: none">-Lane A:0.5 mg Hela Whole Cell Lysate.-Lane B:0.5 mg A549 Whole Cell Lysate.-Lane C:0.5 mg HuT78 Whole Cell Lysate.-2 µL anti-DCUN1D2 rabbit monoclonal antibody and 15 µl of 50 % Protein G agarose.-Primary antibody:-Anti-DCUN1D2 rabbit monoclonal antibody, at 1:100 dilution.-Secondary antibody:-Dylight 800-labeled antibody to rabbit IgG (H+L), at 1:5000 dilution.-Developed using the odyssey technique.-Performed under reducing conditions.-Predicted band size: 30 kDa.-Observed band size: 30 kDa |
| Application: | ELISA,IP,WB |
| Recommended | WB: 1:500-1:2000; ELISA: 1:5000-1:10000; IP: 1-4 µL/mg of lysate |

Properties

| | |
|----------------------|--|
| Stability & Storage: | Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free. |
| Shipping: | Shipping with blue ice. |

Antigen Details

Immunogen: Recombinant Protein: Human DCUN1D2 / DCN1-like protein 2 Protein (TMPY-00580)

Antigen Species: Human

Synonyms: DCN1, defective in cullin neddylation 1, domain containing 2;C13orf17

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481