

Anti-GPD1 Antibody (6D493)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	6D493
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:5000-1:10000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human GPD1 / Glycerolphosphate Dehydrogenase 1 Protein (TMPY-03755)
Antigen Species:	Human
Synonyms:	GPD-C;glycerol-3-phosphate dehydrogenase 1 (soluble);GPDH-C;HTGTI

Research Background

GPD1 (Glycerol-3-Phosphate Dehydrogenase 1) is a Protein Coding gene. 2 alternatively spliced human isoforms have been reported. GPD1 is a member of the NAD-dependent glycerol-3-phosphate dehydrogenase family. The encoded protein plays a critical role in carbohydrate and lipid metabolism by catalyzing the reversible conversion of dihydroxyacetone phosphate (DHAP) and reduced nicotine adenine dinucleotide (NADH) to glycerol-3-phosphate (G3P) and NAD⁺. It also reduces nicotine adenine dinucleotide (NADH) to glycerol-3-phosphate (G3P) and NAD⁺. Meanwhile, GPD1 and mitochondrial glycerol-3-phosphate dehydrogenase also form a glycerol phosphate shuttle that facilitates the transfer of reducing equivalents from the cytosol to mitochondria. Diseases associated with GPD1 include Hypertriglyceridemia, Transient Infantile, and Myopathy, Distal, 1.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481