

Anti-LILRB1/CD85j/ILT2 Antibody-PE (5D841)

Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	PE
Clone:	5D841
Purification:	Protein A

Applications

Verified Activity:	Flow cytometric analysis of Human LILRB1(CD85j) expression on human peripheral blood monocytes. Cells were stained with PE-conjugated Mouse anti-Human LILRB1(CD85j). The histogram were derived from gated events with the forward and side light-scatter characteristics of viable monocytes.
Application:	FCM
Recommended	10 µl/Test, 0.1 mg/ml

Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human CD85 / ILT2 / ILR1 protein (TMPY-03995)
Antigen Species:	Human
Synonyms:	ILT-2;XXbac-BCX85G21.4;LIR-1;ILT2;leukocyte immunoglobulin like receptor B1;CD85j;MIR-7;MIR7;LIR1
Biology Area:	ITIM/ITAM Immunoreceptors and Related Molecules

Research Background

LILRB1 (Leukocyte Immunoglobulin Like Receptor B1) is a Protein Coding gene. The encoded protein belongs to the subfamily B class of LIR receptors which contain two or four extracellular immunoglobulin domains, a transmembrane domain, and two to four cytoplasmic immunoreceptor tyrosine-based inhibitory motifs (ITIMs). LILRB1 is encoded within the leukocyte receptor complex on 19q13.4, previously implicated to be a susceptibility region to systemic lupus erythematosus (SLE). LILRB1 binds to a wide range of Major Histocompatibility Complex Class I molecules (MHCIs) and transduces negative signals that can, for example, prevent the killing of MHC1-expressing cells. LILRB1 transcription is correlated with surface protein expression on NK cells. Diseases associated with LILRB1 include Alzheimer's Disease 12 and Severe Pre-Eclampsia.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481