

## Anti-Growth Hormone Antibody (11177)

### Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	11177
Purification:	Protein A

### Applications

Verified Activity:	Immunochemical staining of human GH1 in human placenta with mouse monoclonal antibody at 1:60 dilution, formalin-fixed paraffin embedded sections.
Application:	IHC-P
Recommended	IHC-P: 1:30-1:100

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human Growth Hormone/GH1 Protein (TMPY-00195)
Antigen Species:	Human
Synonyms:	GH1;IGHD1B;hGH-N;GH;GHB5;GH-N;Growth hormone 1;GHN
Biology Area:	Neurotransmitter Associated Enzymes

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481