

Anti-Growth Hormone Antibody (2A375)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	2A375
Purification:	Protein A

Applications

Verified Activity:	Immunochemical staining of human Growth Hormone in human placenta with rabbit monoclonal antibody at 1:1000 dilution, formalin-fixed paraffin embedded sections.
Application:	IHC-P
Recommended	IHC-P: 1:500-1:3000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Human Growth Hormone / GH1 Protein (TMPY-00195)
Antigen Species:	Human
Synonyms:	GH1;GH-N;IGHD1B;hGH-N;Growth hormone 1;GHB5;GH;GHN
Biology Area:	Neurotransmitter Associated Enzymes

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481