

Anti-Nectin-4 Antibody (60477)

Product Details

| | |
|---------------|--------------|
| Ig Type: | Mouse IgG1 |
| Reactivity: | Human |
| Conjugation: | Unconjugated |
| Clone: | 60477 |
| Purification: | Protein A |

Applications

| | |
|--------------|----------------------|
| Application: | ELISA |
| Recommended | ELISA: 1:1000-1:2000 |

Properties

| | |
|----------------------|--|
| Stability & Storage: | Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------------|--|
| Immunogen: | Recombinant Protein: Human Nectin-4 Protein (TMPY-05537) |
| Antigen Species: | Human |
| Synonyms: | EDSS1;PRR4;LNIR;nectin-4;PVRL4 |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481