

Anti-CD163 Antibody (9L354)

Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	9L354
Purification:	Protein A

Applications

Verified Activity:	Immunochemical staining of human CD163 in human liver with mouse monoclonal antibody at 1:5000 dilution, formalin-fixed paraffin embedded sections.
Application:	IHC-P
Recommended	IHC-P: 1:2000-1:10000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	E. coli-derived Human CD163 fragment
Antigen Species:	Human
Synonyms:	CD163 molecule;sCD163;HbSR;SCAR11;M130;MM130;GHI/61;CD163;RM3/1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481