

Anti-Ebola virus EBOV (subtype Zaire, strain H.sapiens-wt/GIN/2014/Kissidougou-C15)
Nucleoprotein/NP Antibody (6L412)

Product Details

Ig Type: Mouse IgG1
Conjugation: Unconjugated
Clone: 6L412
Purification: Protein A

Applications

Application: ELISA
Recommended: ELISA: 1:1000-1:2000

Properties

Stability & Storage: Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein: Ebola virus EBOV (subtype Zaire, strain H.sapiens-wt/GIN/2014/Kissidougou-C15) Nucleoprotein / NP Protein (TMPY-00312)
Antigen Species: EBOV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481