

Anti-Ebola virus EBOV (Sudan ebolavirus, strain Gulu) Nucleoprotein/NP Antibody (6W419)

Product Details

Ig Type:	Mouse IgG1
Conjugation:	Unconjugated
Clone:	6W419
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:1000-1:2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Ebola virus EBOV (Sudan ebolavirus, strain Gulu) Nucleoprotein / NP Protein (TMPY-00308)
Antigen Species:	EBOV
Synonyms:	Sec61 α 1 subunit (<i>S. cerevisiae</i>); Sec61 alpha 1 subunit (<i>S. cerevisiae</i>)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481