

## Anti-Ebola virus EBOV (subtype Sudan, strain Gulu) VP40 Antibody (6M955)

### Product Details

Ig Type:	Mouse IgG1
Conjugation:	Unconjugated
Clone:	6M955
Purification:	Protein A

### Applications

Verified Activity:	<p>1. Anti-Ebola virus VP40 mouse monoclonal antibody at 1:1000 dilution.</p> <ul style="list-style-type: none"><li>-Sample: Ebola virus EBOV (subtype Sudan, strain Gulu) VP40 Recombinant Protein</li><li>-Lane A: 30ng</li><li>-Lane B: 10ng</li><li>-Secondary</li><li>-Goat Anti-Mouse IgG (H+L)/HRP at 1/10000 dilution.</li><li>-Developed using the ECL technique.</li><li>-Performed under reducing conditions.</li></ul> <p>2. Anti-Ebola virus VP40 mouse monoclonal antibody at 1:1000 dilution</p> <ul style="list-style-type: none"><li>-Sample: Recombinant Protein 30 ng</li><li>-Lane A: Ebola virus EBOV (subtype Zaire, strain H.sapiens-wt/GIN/2014/Kissidougou-C15) Matrix protein VP40 Protein (His Tag)</li><li>-Lane B: Ebola virus EBOV (subtype Sudan, strain Gulu) VP40 / Matrix protein VP40 Protein</li><li>-Lane C: Ebola virus EBOV (subtype Bundibugyo, strain Uganda 2007) VP40 / Matrix protein VP40 Protein</li><li>-Secondary</li><li>-Goat Anti-Mouse IgG (H+L)/HRP at 1/10000 dilution.</li><li>-Developed using the ECL technique.</li><li>-Performed under reducing conditions.</li></ul>
Application:	ELISA,WB
Recommended	WB: 1:1000-1:5000; ELISA: 1:1000-1:2000

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein: EBOV Ebola virus VP40 Protein (TMPY-00216)

Antigen Species: EBOV

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481