

Anti-Ebola virus EBOV (subtype Bundibugyo, strain Uganda 2007) VP40/Matrix protein VP40 Antibody (1H311)

Product Details

Ig Type:	Mouse IgG1
Conjugation:	Unconjugated
Clone:	1H311
Purification:	Protein A

Applications

Verified Activity:	<p>1. Anti-Ebola virus EBOV (subtype Bundibugyo, strain Uganda 2007) VP40 / Matrix protein VP40 mouse monoclonal antibody at 1:1000 dilution.</p> <ul style="list-style-type: none">-Sample:Ebola virus EBOV (subtype Bundibugyo, strain Uganda 2007) VP40 / Matrix protein VP40-Lane A: 30ng-Lane B: 10ng-Secondary-Goat Anti-Mouse IgG H&L (Dylight800) at 1/15000 dilution.-Developed using the Odyssey technique.-Performed under reducing conditions.
	<p>2. Anti-Ebola virus EBOV (subtype Bundibugyo, strain Uganda 2007) VP40 / Matrix protein VP40 mouse monoclonal antibody at 1:1000 dilution.</p> <ul style="list-style-type: none">-Sample: Recombinant Protein 10 ng-Lane A: Ebola virus EBOV (subtype Zaire, strain H.sapiens-wt/GIN/2014/Kissidougou-C15) Matrix protein VP40 Protein-Lane B: Ebola virus EBOV (subtype Sudan, strain Gulu) VP40 / Matrix protein VP40 Protein-Lane C: Ebola virus EBOV (subtype Bundibugyo, strain Uganda 2007) VP40 / Matrix protein VP40 Protein-Secondary-Goat Anti-Mouse IgG (H+L)/HRP at 1/10000 dilution.-Developed using the ECL technique.-Performed under reducing conditions.
Application:	ELISA,WB
Recommended	WB: 1:1000-1:5000; ELISA: 1:1000-1:2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein: Ebola virus EBOV (subtype Bundibugyo, strain Uganda 2007) VP40 / Matrix protein VP40 Protein (TMPY-00065)
Antigen Species: EBOV
Synonyms: Matrix protein VP40

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481