

Anti-Dengue virus (DENV) (type 2, strain New Guinea C/PUO-218 hybrid) E/Envelope Protein (Domain III) Antibody (5L714)

Product Details

Ig Type:	Mouse IgG1
Conjugation:	Unconjugated
Clone:	5L714
Purification:	Protein A

Applications

Verified Activity:	Anti-Dengue virus (DENV) (type 2, strain New Guinea C/PUO-218 hybrid) E / Envelope Protein (Domain III) rabbit polyclonal antibody at 1:2000 dilution -Sample:Dengue virus (DENV) (type 2, strain New Guinea C/PUO-218 hybrid) E / Envelope Protein (Domain III) Recombinant Protein -Lane A: 30ng -Secondary -Rabbit Anti-Mouse IgG (H+L)/HRP at 1/10000 dilution -Performed under reducing conditions.
Application:	WB
Recommended	WB: 1:1000-1:5000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Dengue virus (DENV) (type 2, strain New Guinea C/PUO-218 hybrid) E / Envelope Protein (Domain III) Protein (TMPY-04700)
Antigen Species:	DENV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481