

Anti-Zika virus (ZIKV) (strain Zika SPH2015) ZIKV-NS1 Antibody (4I627)

Product Details

Ig Type:	Mouse IgG1
Conjugation:	Unconjugated
Clone:	4I627
Purification:	Protein A

Applications

Application:	ELISA(Det)
--------------	------------

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Zika virus (ZIKV) (strain Zika SPH2015) ZIKV-NS1 Protein (TMPY-04870)
Antigen Species:	ZIKV

Research Background

Zika virus NS1 antigen is one of seven non-structural proteins. NS1 is involved in RNA replication. The possible effects of NS1 on hosts include: localizes to host cell surface and secreted extracellularly, modulates signalling of the innate immune system, has possible damages to platelets and endothelial cells through anti-NS1 antibodies.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481