

Anti-Vesicular stomatitis Indiana virus (VSIV) (strain San Juan) glycoprotein/GP Antibody (6S617)

Product Details

Ig Type: Rabbit IgG
Conjugation: Unconjugated
Clone: 6S617
Purification: Protein A

Applications

Application: ELISA
Recommended: ELISA: 1:5000-1:10000

Properties

Stability & Storage: Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein: Vesicular stomatitis Indiana virus (VSIV) (strain San Juan) glycoprotein / GP Protein
Antigen Species: VSIV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481