

Anti-Monkeypox Virus (MPXV) A29 Antibody (9W603)

Product Details

Ig Type:	Mouse IgG1
Conjugation:	Unconjugated
Clone:	9W603
Purification:	Protein A

Applications

Verified Activity:	Anti-Monkeypox Virus (MPXV) Protein A29 mouse monoclonal antibody at 1:1000 dilution. -Lane A: Monkeypox Virus (MPXV) Protein A29 10ng -Lane B: Vaccinia virus (VACV) (strain Copenhagen) A27L Protein 10ng -Lane C: Monkeypox Virus (MPXV) (Clade Ib) Protein A29L (10ng) -Secondary -Goat Anti-Mouse IgG (H+L)/HRP at 1/10000 dilution. -Developed using the ECL technique. -Performed under reducing conditions.
Application:	ELISA,WB
Recommended	WB: 1:1000-1:5000; ELISA: 1:1000-1:2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Monkeypox Virus (MPXV) Protein A29 (TMPY-06522)
Antigen Species:	MPXV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481