

## Anti-Varicella-zoster virus (VZV)/HHV-3 Glycoprotein E/gE Antibody (2E390)

## Product Details

Ig Type:	Mouse IgG1
Conjugation:	Unconjugated
Clone:	2E390
Purification:	Protein A

## Applications

Verified Activity:	Anti-VZV / HHV-3 gE / glycoprotein E mouse monoclonal antibody in 1 µg/mL. -Sample: Recombinant Protein 10 ng -Lane A: VZV / HHV-3 gE / glycoprotein E Protein -Lane B: VZV (strain Oka vaccine) Glycoprotein B Protein -Lane C: VZV (strain Oka vaccine) gH & gL Heterodimer Protein -Secondary -Goat Anti-Mouse IgG (H+L)/HRP at 1/10000 dilution. -Developed using the ECL technique. -Performed under reducing conditions.
Application:	ELISA,WB
Recommended	WB: 0.2-1 µg/mL; ELISA: 0.5-1 µg/mL

## Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	Recombinant Protein: Varicella-zoster virus (VZV) / HHV-3 gE / glycoprotein E Protein (TMPY-06357)
Antigen Species:	VZV

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481