

Anti-Varicella-zoster virus (VZV)/HHV-3 Glycoprotein E/gE Antibody (3V671)

Product Details

Ig Type:	Rabbit IgG
Conjugation:	Unconjugated
Clone:	3V671
Purification:	Protein A

Applications

Verified Activity:	Anti-VZV / HHV-3 gE / glycoprotein E rabbit monoclonal antibody at 0.5 µg/mL. -Sample: Recombinant Protein 5 ng -Lane A: VZV / HHV-3 gE / glycoprotein E Protein -Lane B: VZV (strain Oka vaccine) gB Protein -Lane C: VZV (strain Oka vaccine) gH & gL Protein -Secondary -Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution. -Developed using the ECL technique. -Performed under reducing conditions.
Application:	ELISA,WB
Recommended	WB: 0.1-1 µg/mL; ELISA: 0.1-0.2 µg/mL

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Varicella-zoster virus (VZV) / HHV-3 gE / glycoprotein E Protein (TMPY-06357)
Antigen Species:	VZV

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481