

## Anti-FCGR4 Antibody-APC (7C860)

### Product Details

Ig Type:	Rabbit IgG
Reactivity:	Mouse
Conjugation:	APC
Clone:	7C860
Purification:	Protein A

### Applications

Verified Activity:	Flow cytometric analysis of FCGR4 expression on Raw264.7 cells.
Application:	FCM
Recommended	10 µl/Test, 0.1 mg/ml

### Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight. Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Mouse FCGR4 protein (TMPY-01087)
Antigen Species:	Mouse
Synonyms:	FcyRIV;Fcrl3;Fcgr3a;FcγRIV;CD16-2;4833442P21Rik;FcγammaRIV;Fc receptor, IgG, low affinity IV
Biology Area:	Fc Receptors

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481