

Anti-FCGR4 Antibody-PerCP (7C860)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Mouse
Conjugation:	PerCP
Clone:	7C860
Purification:	Protein A

Applications

Verified Activity:	Flow cytometric analysis of FCGR4 expression on Raw264.7 cells.
Application:	FCM
Recommended	10 µl/Test, 0.1 mg/ml

Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight. Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Mouse FCGR4 protein (TMPY-01087)
Antigen Species:	Mouse
Synonyms:	FcgRIV;Fcgr3a;4833442P21Rik;FcγRIV;Fcrl3;Fc receptor, IgG, low affinity IV;CD16-2; FcγgammaRIV
Biology Area:	Fc Receptors

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481