

Anti-Meprin alpha/MEP1A Antibody (2Q198)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Mouse
Conjugation:	Unconjugated
Clone:	2Q198
Purification:	Protein A

Applications

Verified Activity:	1. Immunochemical staining of mouse MEP1A in mouse kidney with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections). 2. Immunochemical staining of mouse MEP1A in mouse intestine with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections). The left panel: tissue incubated with primary antibody; The right panel: tissue incubated with the mixture of primary antibody and antigen (recombinant protein).
Application:	ELISA,IHC-P
Recommended	ELISA: 1:5000-1:10000; IHC-P: 1:100-1:500

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Mouse MEP1A Protein (TMPY-01318)
Antigen Species:	Mouse
Synonyms:	meprin A, α (PABA peptide hydrolase);PPHA;Meprin α /MEP1A;meprin A, alpha (PABA peptide hydrolase)

Research Background

Meprin A subunit alpha, also known as MEP1A, and Endopeptidase-2, is a single-pass type I membrane protein that belongs to the peptidase M12A family. MEP1A contains one EGF-like domain, one MAM domain, and one MATH domain. Meprins are unique plasma membrane and secreted metalloproteinases that are highly regulated at the transcriptional and post-translational levels. Meprin alpha and beta subunits are abundantly expressed in kidney and intestinal epithelial cells, are secreted into the urinary tract and intestinal lumen and are found in leukocytes and cancer cells under certain conditions. Meprins are capable of proteolytically degrading extracellular matrix proteins, processing bioactive proteins, and play a role in inflammatory processes. Meprin A and B are highly regulated, secreted and cell-surface homo- and hetero-oligomeric enzymes. Meprins are abundantly expressed in the kidney and intestine. The multidomain alpha and beta subunits have high sequence identity. They have very different substrate specificities, oligomerization potentials, and are differentially regulated. Meprin A appears to be an important therapeutic target and urinary excretion appears to be a potential biomarker of acute kidney injury (AKI).

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481