

Anti-SLAMF6 Antibody (6R165)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Mouse
Conjugation:	Unconjugated
Clone:	6R165
Purification:	Protein A

Applications

	Anti-SLAMF6 rabbit monoclonal antibody at 1:500 dilution. -Lane A: Mouse heart tissue lysate. -Lysates/proteins at 30 µg per lane. -Secondary
Verified Activity:	-Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution. -Developed using the ECL technique. -Performed under reducing conditions. -Predicted band size:39 kDa. -Observed band size:44 kDa
Application:	WB
Recommended	WB: 1:500-1:2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Mouse Ly108/SLAMF6 Protein (TMPY-01809)
Antigen Species:	Mouse
Synonyms:	SLAM family member 6;Ly108;KALI;SF2000;NTB-A;KALIB;CD352;NTBA
Biology Area:	ITIM/ITAM Immunoreceptors and Related Molecules

Research Background

SLAM family member 6, also known as Activating NK receptor, NK-T-B-antigen, NTB-A, SLAMF6, KALI and Ly18, is a single-pass type I membrane protein that belongs to the CD2 subfamily of the immunoglobulin superfamily. SLAMF6 / Ly18 contains one Ig-like (immunoglobulin-like) domain. It is expressed by all (resting and activated) natural killer cells (NK), T- and B-lymphocytes. SLAMF6 / Ly18 triggers cytolytic activity only in natural killer cells (NK) expressing high surface densities of natural cytotoxicity receptors. SLAMF6 / Ly18 is a homodimer. It interacts with PTN6 and, upon phosphorylation, with PTN11 and SH2D1A/SAP. SLAMF6 / Ly18 undergoes tyrosine phosphorylation and associates with the Src homology 2 domain-containing protein (SH2D1A) as well as with SH2 domain-containing phosphatases (SHPs). It may function as a coreceptor in the process of NK cell activation. SLAMF6 / Ly18 can also mediate inhibitory signals in NK cells from X-linked lymphoproliferative patients.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481