

## Anti-Serpin A3c Antibody (9E761)

### Product Details

Ig Type:	Rabbit IgG
Reactivity:	Mouse
Conjugation:	Unconjugated
Clone:	9E761
Purification:	Protein A

### Applications

Application:	ELISA,WB
Recommended	ELISA: 1:5000-1:10000; WB: 1:500-1:2000

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Mouse SerpinA3c Protein (TMPY-01329)
Antigen Species:	Mouse
Synonyms:	Serpina3c;1A1;serine (or cysteine) peptidase inhibitor, clade A, member 3C;spi2;Klkbp;Kalbp; Serpin A3C;MGC144822

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481