

## Anti-CD229 Antibody-PerCP (7Q20)

### Product Details

Ig Type:	Rabbit IgG
Reactivity:	Mouse
Conjugation:	PerCP
Clone:	7Q20
Purification:	Protein A

### Applications

Verified Activity:	Analysis of Ly9 (CD229) expression on spleen lymphocytes.
Application:	FCM
Recommended	5 µl/Test, 0.1 mg/ml

### Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight. Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Mouse Ly9 / CD229 / SLAMF3 protein (TMPY-01804)
Antigen Species:	Mouse
Synonyms:	T100;SLAMF3;lymphocyte antigen 9;Ly-9;AI893573;CD229;Lgp100
Biology Area:	ITIM/ITAM Immunoreceptors and Related Molecules

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481