

## Anti-CD229 Antibody-PE (7Q20)

### Product Details

|               |            |
|---------------|------------|
| Ig Type:      | Rabbit IgG |
| Reactivity:   | Mouse      |
| Conjugation:  | PE         |
| Clone:        | 7Q20       |
| Purification: | Protein A  |

### Applications

|                    |   |
|--------------------|---|
| Verified Activity: | Analysis of Ly9 (CD229) expression on spleen lymphocytes. |
| Application:       | FCM   |
| Recommended        | 2 µl/Test, 0.1 mg/ml                                      |

### Properties

|                      |  |
|----------------------|--|
| Stability & Storage: | Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight. Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal. |
| Shipping:            | Shipping with blue ice.  |

### Antigen Details

|                  |  |
|------------------|--|
| Immunogen:       | Recombinant Protein: Mouse Ly9 / CD229 / SLAMF3 protein (TMPY-01804) |
| Antigen Species: | Mouse  |
| Synonyms:        | SLAMF3;Ly-9;lymphocyte antigen 9;Lgp100;T100;AI893573;CD229          |
| Biology Area:    | ITIM/ITAM Immunoreceptors and Related Molecules                      |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481