

## Anti-Interferon alpha 14/IFNA14 Antibody (4Q295)

### Product Details

|               |              |
|---------------|--------------|
| Ig Type:      | Rabbit IgG   |
| Reactivity:   | Mouse        |
| Conjugation:  | Unconjugated |
| Clone:        | 4Q295        |
| Purification: | Protein A    |

### Applications

|              |                       |
|--------------|-----------------------|
| Application: | ELISA                 |
| Recommended  | ELISA: 1:5000-1:10000 |

### Properties

|                      |  |
|----------------------|--|
| Stability & Storage: | Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free. |
| Shipping:            | Shipping with blue ice.  |

### Antigen Details

|                  |   |
|------------------|---|
| Immunogen:       | Recombinant Protein: Mouse IFNA14 Protein (TMPY-03815)          |
| Antigen Species: | Mouse   |
| Synonyms:        | interferon, alpha 14;interferon, $\alpha$ 14;OTTMUSG00000007662 |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481