

## Anti-ST6GALNAC2 Antibody (2T197)

### Product Details

Ig Type:	Rabbit IgG
Reactivity:	Mouse
Conjugation:	Unconjugated
Clone:	2T197
Purification:	Protein A

### Applications

Application:	ELISA
Recommended	ELISA: 1:5000-1:10000

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Mouse ST6GALNAC2 protein (TMPY-02685)
Antigen Species:	Mouse
Synonyms:	ST6-N-acetylgalactosaminide $\alpha$ -2,6-sialyltransferase 2;ST6GalNAc;Siat7;Siat7b;ST6-N-acetylgalactosaminide alpha-2,6-sialyltransferase 2;II

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481