

Anti-CREG1 Antibody (20152)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Mouse
Conjugation:	Unconjugated
Clone:	20152
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:5000-1:10000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Mouse CREG1 Protein (TMPY-03420)
Antigen Species:	Mouse
Synonyms:	CREG;cellular repressor of E1A-stimulated genes 1

Research Background

CREG1 belongs to the CREG family. It is a adenovirus E1A protein which both activates and represses gene expression to promote cellular proliferation and inhibit differentiation. Thus it may contribute to the transcriptional control of cell growth and differentiation. The transcriptional control activity of cell growth requires interaction with IGF2R. CREG1 also antagonizes transcriptional activation and cellular transformation by E1A. It shares limited sequence similarity with E1A and binds both the general transcription factor TBP and the tumor suppressor pRb in vitro. CREG1 gene may contribute to the transcriptional control of cell growth and differentiation.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481