

Anti-IL-28B Antibody (5K886)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Mouse
Conjugation:	Unconjugated
Clone:	5K886
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:5000-1:10000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Mouse IL-28B / IFN-lambda-3 Protein (TMPY-00454)
Antigen Species:	Mouse
Synonyms:	IL28C;interferon, λ3;IL-28B;IL28B;interferon, lambda 3

Research Background

Interleukin-28B (IL-28B) also known as Interferon lambda-3 and IFN-lambda-3, belongs to the type III interferon family of cytokines and are highly similar to IL-29. IL-28B belongs to the newly described interferon lambda (IFNλ) family of cytokines. IL-28B is a cytokine with immunomodulatory activity. It functions in Up-regulating MHC class I antigen expression. IL-28B displays potent antiviral activity and antitumor activity. This cytokine serves as a ligand for the heterodimeric class II cytokine receptor composed of IL10RB and IL28RA. The ligand/receptor complex seems to signal through the Jak-STAT pathway. IL-28B, like IL-12, is capable of robustly enhancing adaptive immunity. Moreover, we describe for the first time how IL-28B reduces regulatory T-cell populations during DNA vaccination, whereas IL-12 increases this cellular subset. We also show that IL-28B, unlike IL-12, can increase the percentage of splenic CD8+ T cells in vaccinated animals and that these cells are more granular and have higher antigen-specific cytolytic degranulation compared with cells taken from animals that received IL-12 as an adjuvant.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481