

Anti-Albumin Antibody (7N226)

Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	7N226
Purification:	Protein A

Applications

Verified Activity:	Anti-Albumin rabbit polyclonal antibody at 1:500 dilution. -Lane A: HepG2 Whole Cell lysate. -Lysates/proteins at 30 µg per lane. -Secondary -Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution. -Developed using the ECL technique. -Performed under reducing conditions. -Predicted band size:69 kDa
Application:	ELISA,WB
Recommended	WB: 1:500-1:2000; ELISA: 1:1000-1:2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Human Albumin
Antigen Species:	Human
Synonyms:	PRO1341;albumin;PRO0883;ANALBA;FDAH;PRO0903

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481