

Anti-Human Serum Albumin/HSA/ALB Antibody-FITC (40456)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	FITC
Clone:	40456
Purification:	Protein A

Applications

Verified Activity:	Flow cytometric analysis of Human HSA expression on HepG2 cells. The cells were treated according to manufacturer's manual (BD Pharmingen™ Cat. No. 554714), and then stained with FITC-conjugated anti-Human HSA. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.
Application:	FCM
Recommended	10 µl/Test, 0.1 mg/ml

Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Human Serum Albumin / HSA
Antigen Species:	Human
Synonyms:	PRO1341;PRO2675;PRO0883;PRO2619;PRO1708;FDAHT;GIG42;GIG20;PRO0903;UNQ696/PRO1341;PRO2044;HSA

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481