

## Anti-Human Serum Albumin/HSA/ALB Antibody-HRP (40456)

### Product Details

Ig Type:	Rabbit IgG
Reactivity:	Human
Conjugation:	HRP
Clone:	40456
Purification:	Protein A

### Applications

Application:	ELISA
Recommended	ELISA: 0.1-1 µg/mL

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Human Serum Albumin / HSA
Antigen Species:	Human
Synonyms:	PRO0883;GIG42;UNQ696/PRO1341;PRO2675;PRO1341;PRO2044;FDAHT;GIG20;HSA;PRO2619;PRO0903;PRO1708

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481