

Anti-CD45RA Antibody-FITC (6E481)

Product Details

| | |
|--------------|-------------|
| Ig Type: | Mouse IgG2a |
| Reactivity: | Human |
| Conjugation: | FITC |
| Clone: | 6E481 |

Applications

| | |
|--------------------|--|
| Verified Activity: | Flow cytometric analysis of Human CD45RA expression on human whole blood lymphocytes. Cells were stained with FITC-conjugated anti-Human CD45RA. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable lymphocytes. |
| Application: | FCM |
| Recommended | 2 µl/Test, 0.1 mg/ml |

Properties

| | |
|----------------------|--|
| Stability & Storage: | Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight. Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------------|---|
| Immunogen: | 0 |
| Antigen Species: | Human |
| Synonyms: | protein tyrosine phosphatase, receptor type, C;CD45 antigen;B220;receptor-type tyrosine-protein phosphatase C;L-CA;CD45;LY5;LCA;PTPRC;CD45R |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481