

## Anti-CD45RA Antibody-PE (6E481)

### Product Details

Ig Type:	Mouse IgG2a
Reactivity:	Human
Conjugation:	PE
Clone:	6E481

### Applications

Verified Activity:	Flow cytometric analysis of Human CD45RA expression on human whole blood lymphocytes. Cells were stained with PE-conjugated anti-Human CD45RA. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable lymphocytes.
Application:	FCM
Recommended	0.5 µl/Test, 0.1 mg/ml

### Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight. Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	0
Antigen Species:	Human
Synonyms:	LY5;CD45;CD45R;B220;L-CA;PTPRC;CD45 antigen;receptor-type tyrosine-protein phosphatase C;LCA;protein tyrosine phosphatase, receptor type, C

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481