

Anti-CD45RA Antibody (1M474)

Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	1M474
Purification:	Protein A

Applications

Verified Activity:	Flow cytometric analysis of Human CD45RA expression on human whole blood lymphocytes. Cells were stained with purified anti-Human CD45RA, then a FITC-conjugated second step antibody. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable lymphocytes.
Application:	FCM
Recommended	FCM: 1:25-1:100

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	0
Antigen Species:	Human
Synonyms:	PTPRC;B220;LCA;L-CA;CD45R;receptor-type tyrosine-protein phosphatase C;LY5;CD45 antigen; CD45;protein tyrosine phosphatase, receptor type, C

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481