

## Anti-CD45RA Antibody-PerCP (1M474)

### Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	PerCP
Clone:	1M474
Purification:	Protein A

### Applications

Verified Activity:	Flow cytometric analysis of Human CD45RA expression on human whole blood lymphocytes. Cells were stained with PerCP-conjugated anti-Human CD45RA. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable lymphocytes.
Application:	FCM
Recommended	10 µl/Test, 0.1 mg/ml

### Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	0
Antigen Species:	Human
Synonyms:	B220;CD45;L-CA;protein tyrosine phosphatase, receptor type, C;PTPRC;CD45 antigen;CD45R;LY5;LCA;receptor-type tyrosine-protein phosphatase C

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481