

Anti-CD45RA Antibody-FITC (1M474)

Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	FITC
Clone:	1M474
Purification:	Protein A

Applications

Verified Activity:	Flow cytometric analysis of Human CD45RA expression on human whole blood lymphocytes. Cells were stained with FITC-conjugated anti-Human CD45RA. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable lymphocytes.
Application:	FCM
Recommended	5 µl/Test, 0.1 mg/ml

Properties

Stability & Storage:	Store at 2°C-8°C for 12 months, do not freeze. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	0
Antigen Species:	Human
Synonyms:	LCA;PTPRC;CD45 antigen;protein tyrosine phosphatase, receptor type, C;LY5;receptor-type tyrosine-protein phosphatase C;CD45R;B220;L-CA;CD45

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481