

## Anti-SARS-CoV-2 Nucleocapsid Antibody (6H645)

## Product Details

|               |              |
|---------------|--------------|
| Ig Type:      | Mouse IgG1   |
| Conjugation:  | Unconjugated |
| Clone:        | 6H645        |
| Purification: | Protein A    |

## Applications

Anti-SARS-CoV-2 Nucleocapsid mouse monoclonal antibody at 1:1000 dilution.

-Sample: SARS-CoV-2 (2019-nCoV) Nucleocapsid-His recombinant Protein

-Lane A: 10ng

Verified Activity:

-Secondary

-Goat Anti-Mouse IgG (H+L)/HRP at 1/10000 dilution.

-Developed using the ECL technique.

-Performed under reducing conditions.

Application: ELISA,WB

Recommended WB: 1:500-1:2000; ELISA: 1:1000-1:2000

## Properties

Stability & Storage: Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.

Shipping: Shipping with blue ice.

## Antigen Details

Immunogen: Recombinant Protein: SARS-CoV-2 Nucleocapsid CTD Protein (TMPY-06113)

Antigen Species: SARS-CoV-2

## Research Background

Coronaviruses are enveloped viruses with a positive-sense RNA genome and with a nucleocapsid of helical symmetry. Coronavirus nucleoproteins localize to the cytoplasm and the nucleolus, a subnuclear structure, in both virus-infected primary cells and in cells transfected with plasmids that express N protein. The coronavirus N protein is required for coronavirus RNA synthesis and has RNA chaperone activity that may be involved in template switch. Nucleocapsid protein is the most abundant protein of coronavirus. During virion assembly, N protein binds to viral RNA and leads to the formation of the helical nucleocapsid. Nucleocapsid protein is a highly immunogenic phosphoprotein also implicated in viral genome replication and in modulating cell signaling pathways. Because of the conservation of the N protein sequence and its strong immunogenicity, the N protein of coronavirus is chosen as a diagnostic tool.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481