

## Anti-MMAE Antibody (20521)

## Product Details

Ig Type:	Mouse IgG1
Reactivity:	other
Conjugation:	Unconjugated
Clone:	20521
Purification:	Protein A

## Applications

Verified Activity:	<p>1. Immobilized MMAE-ADC, Naked Antibody, MMAE-OVA, MMAF-OVA, OVA, at 1 µg/mL (100 µL/well), then add the Anti-MMAE Antibody to detect the cross-reactivity. (Validation Experiment)</p> <p>2. Immobilized Anti-MMAE Antibody at 2 µg/mL (100 µL/well) can bind the MMAE-ADC, the EC50 is 50-1000 ng/mL. The naked antibody (without MMAE) is the negative control.</p> <p>3. Immobilized Anti-MMAE Antibody at 2 µg/mL (100 µL/well), gradient-diluted free MMAE from 10 µg/ml (50 µL/well) blocked the binding of MMAE Antibody to MMAE-HRP 2 µg/mL (50 µL/well). (Validation Experiment)</p>
Application:	ELISA
Recommended	ELISA: 1:1000-1:2000

## Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	MMAE conjugated to BSA
Antigen Species:	other

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481