

## Anti-FAM19A5 Antibody (8Z48)

### Product Details

Ig Type:	Mouse IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	8Z48
Purification:	Protein A

### Applications

Application:	ELISA
Recommended	ELISA: 1:1000-1:2000

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Human TFAA5/FAM19A5 Protein
Antigen Species:	Human
Synonyms:	Tafa5;Fam19a5;UNQ5208;Chemokine-like protein TFAA-5;QLLK5208

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481