

Anti-GHITM/MICS1 Antibody (2V609)

Product Details

Ig Type:	Mouse IgG2a
Reactivity:	Dermatophagoides pteronyssinus
Conjugation:	Unconjugated
Clone:	2V609
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:1000-1:2000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Dermatophagoides pteronyssinus GHITM / MICS1 Protein (TMPY-04945)
Antigen Species:	Dermatophagoides pteronyssinus

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481