

## Anti-Der p 1 Protein, Dermatophagoides pteronyssinus Antibody (7K321)

### Product Details

Ig Type:	Rabbit IgG
Reactivity:	Dermatophagoides pteronyssinus
Conjugation:	Unconjugated
Clone:	7K321
Purification:	Protein A

### Applications

Application:	ELISA
Recommended	ELISA: 1:5000-1:10000

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: Dermatophagoides pteronyssinus Der p 1 Protein (TMPY-04960)
Antigen Species:	Dermatophagoides pteronyssinus

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481