

Anti-SOST Antibody (4J191)

Product Details

Ig Type:	Rabbit IgG
Reactivity:	Rat
Conjugation:	Unconjugated
Clone:	4J191
Purification:	Protein A

Applications

Application:	ELISA
Recommended	ELISA: 1:5000-1:10000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Preservative-Free.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: Rat Sclerostin/SOST Protein (TMPY-01633)
Antigen Species:	Rat
Synonyms:	DAND6;CDD;SOST1;VBCH;sclerostin

Research Background

Sclerostin, the protein product of the SOST gene, is a potent inhibitor of bone formation. Sclerostin protein is widely expressed at low levels with highest levels in bone, cartilage, kidney, liver, bone marrow and primary osteoblasts differentiated for 21 days, and was originally identified as an important regulator of bone remodeling, homeostasis, and links bone resorption and bone apposition. Recent studies have revealed that Sclerostin protein inhibits the bone growth probably by binding to the extracellular domain of the Wnt coreceptors LRP5 and LRP6 and disrupting Wnt-induced Frizzled-LRP complex formation.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481